

**PCR genotyping protocol for Rosa26^{-228-DR5-TA-Cerulean}
and Rosa26^{-228-DR5-TA-Cerulean} alleles**

PCR reagents:

1. oligo primers at 20 μ M
Amplification across lox2272 site: Rosa26.S11 + Rosa26.S2 yield a 523 bp Rosa26^{wt} amplicon & a 588 bp Rosa26^{-228-DR5-TA-Cer(-H)} amplicon
Rosa26.S11 5'- CGTGCTGAGCCAGACCTCCAT -3' (bottom)
Rosa26.S2 5'- TCACAAGCAATAATAACCTGTAGT -3' (bottom)
2. Perkin Elmer PCR buffer with MgCl₂
3. 1.25mM dNTP premix (dNTP premix is made using 100 mM NEB dNTP's. The premix contains 250 μ l of each dNTP - A, C, G, & T - and 19 ml sterile water and is stable at -20°C. I freeze this in 1 ml aliquots and thaw and refreeze as needed.)
4. Perkin Elmer Amplitaq Gold
5. genomic DNA samples diluted to 50 ng/ μ l with sterile water

PCR reaction mixture:

15.8 μ l sterile water
2.5 μ l 10X PCR buffer
4 μ l dNTP premix
0.75 μ l primer #1 (Rosa26.S1 or Rosa26.S10)
0.75 μ l primer #2 (Rosa26.S11 or Rosa26.S2)
1 μ l dil. DNA template
0.2 μ l Amplitaq Gold
25 μ l total volume

Cycling conditions:

1 cycle - 94°C x 6 min.
40 cycles - 94°C x 1 min., 60°C x 30 sec., 72°C x 45 sec.
1 cycle - 72°C x 7 min.
hold at 4°C.

Analysis of PCR products:

Load 10 μ l aliquots of reactions + 2 μ l gel loading buffer in a 1% mini-agarose gel (using a 15-well comb). Run gel approximately half-way down.

