


## GENOTYPING PROTOCOL: Myo5B G519R

Investigator: James Goldenring Genome Edit: *Myo5B*<sup>G519R</sup> Allele name: *Myo5B*<sup>em1Jrgo</sup>

## Overview of genotyping strategy:



## Assay for detection of Myo5B G519R allele:

This assay will not differentiate between heterozygous and homozygous animals. Nucleotides in red are specific to the designed Myo5B G519R edited allele.

PCR Primers: Myo5BG519RFwd: CCTAAGCGGACCGACCAGAAC

Myo5BG519RRev1: GTGCCTCGTCTTCTCAACAGTCTG

PCR product = 307 bp

Components	50 ul reaction	Final concentration	PCR program*
5x Phusion Reaction Buffer (NEB #M0530S)	10.0 μL	1X	98°C, 30 seconds
10 mM dNTPs	1.0 µL	200 μΜ	98°C, 10 seconds
10 μM Primer forward	2.5 µL	0.5 μΜ	67°C, 10 seconds
10 μM Primer reverse	2.5 µL	0.5 μΜ	72°C, 30 seconds
Phusion DNA Polymerase (NEB #M0530S)	0.5 µL	0.02 U/μl	Go to 2, 35 X
Nuclease-free water	32.5 µL		72°C, 2 minutes
Genomic DNA	1.0 µL	Less than 1 µg	4°C, ∞



## Assay for detection of Myo5B G519R allele:

This PCR assay amplifies both WT and Myo5B G519R alleles. Digestion with Bpu10I can be used to detect heterozygous versus homozygous Myo5B G519R animals.

PCR Primers: Myo5BG519RFwd1: GAGCGTCTTCAGACAACTCCACA

Myo5BG519RRev1: GTGCCTCGTCTTCTCAACAGTCTG

PCR product = 614 bp

Components	50 ul reaction	Final concentration	PCR program*
5x Phusion Reaction Buffer (NEB #M0530S)	10.0 µL	1X	98°C, 30 seconds
10 mM dNTPs	1.0 µL	200 μΜ	98°C, 10 seconds
10 μM Primer forward	2.5 μL	0.5 μΜ	67°C, 10 seconds
10 μM Primer reverse	2.5 μL	0.5 μΜ	72°C, 30 seconds
Phusion DNA Polymerase (NEB #M0530S)	0.5 μL	0.02 U/μl	Go to 2, 35 X
Nuclease-free water	32.5 µL		72°C, 2 minutes
Genomic DNA	1.0 µL	Less than 1 µg	4°C, ∞

If you use Phusion HF polymerase you can add 0.5 ul Bpu10l enzyme directly to 10  $\mu$ l PCR product and incubate for 1 hour at 37°C. Use of other polymerases may require purification of the PCR product prior to digestion using the Bpu10l buffer provided by NEB. Run digested PCR samples on a 2% or greater agarose gel to separate digest products. Optimization may be required to reliably assess zygosity.

