

TetO-EGFR Genotyping

Primer name	Sequence (5' to 3')
EGFR3421	ACTGTCCAGCCCACCTGTGT
mp-1R	GCCTGCGACGGCGGCATCTGC

PCR Amplification Protocol

1. Denaturation

- Heat the reaction mixture to **95°C** for **5 minutes** to denature the DNA.

2. Cycling Phase (Repeat for **35 cycles**):

- Denaturation:** Heat to **95°C** for **30 seconds**.
- Annealing:** Lower the temperature to **58°C** for **30 seconds** to allow primers to bind.
- Extension:** Increase the temperature to **72°C** for **30 seconds** to allow DNA polymerase to synthesize new strands.

3. Final Extension

- Maintain the reaction at **72°C** for **5 minutes** to complete strand synthesis.

Transgenic Confirmation by Southern Blotting

- Extract **tail DNA** from transgenic founder specimens.
- Digest the DNA using **XhoI-SalI** restriction enzymes.
- Perform **Southern blotting** using a **~500 bp XhoI-SalI probe** targeting the **5' end of the construct** to confirm transgene integration.