

Grm7-I154T Mouse Genotyping

1. Isolate DNA using NaOH method
 - a. Add 300ul of 50mM NaOH to tail and heat at 95C with occasional vortexing until tail is dissolved
 - b. Add 30ul Tris HCl
2. Set up PCR as below:
 - a. Include water (blank) and HET (positive) controls

Reagent	1 Reaction (ul)	MasterMix
5X HF Phusion Buffer Vortex well before use!	4	
5mM dNTPs	1	
10uM Primer 12	1	
10uM Primer 15	1	
DNA	2	n/a
Water	10.8	
Phusion HF polymerase	0.2	
Total Volume	20	

3. Run PCR program:

Temp	Time		<u>Primer sequences:</u> NF12: TCCTGGACACTTGTTCCAGGGACAC NF15: CATGAAGTCCAAACCAGCTTTT <u>Product size:</u> 285
98	0:30		
98	0:10	35 cycles	
60	0:20		
72	0:30		
72	7:00		
4	Hold		

4. Use Zymo kit to purify PCR product- elute in 15 ul of H₂O
5. Set up a BsrI digest:

Reagent	1 Reaction (ul)	MasterMix
PCR product	10	n/a
BsrI	1	
Buffer 3.1	2	
Water	7	
Total Volume	20	

6. Incubate in PCR machine at 65C for 15 minutes.
7. Add DNA loading dye and run 15ul of product on a 3% DNA gel (small comb) for 15-20 min at 150V
8. Expected results:
 - a. WT (+/+): 285 bp
 - b. HET (+/I154T): 285, 150, 135 bp
 - c. HOMO (I154T/I154T): 150, 135 bp

