

PCR genotyping Protocol for *Pdx1^{Cre}* mouse lines.

PCR reagents:

1. oligo primers at 20 μ M

Amplification across lox71/66 (5') site:

Pdx1.5b 5'- TGA GAT TGT ATA TTG CGG TGC A -3'

Pdx1.S3 5'- AGA ACT CTC TAG ATA GTA AGT -3'

Pdx1.5b + Pdx1.S3 yields a **573 bp TM** (targeted mutant allele) band.

Amplification across lox2272 site (3'):

Pdx1.S5 5'- TGA GCA ATT CCA AGC AGC TGG A -3'

Pdx1.3c 5'- ACC TTG CAG TCC TTC TGA AGT -3'

Pdx1.S5 + Pdx1.3c yields a **426bp WT** (wild type allele) band and a **481bp TM** (targeted mutant allele) band.

2. Perkin Elmer PCR buffer with MgCl₂
3. 1.25mM dNTP premix (dNTP premix is made by using 100 mM NEB dNTP's. The premix contains 250 μ l of each dNTP - A, C, G, & T - and 19 ml sterile water and is stable at -20°C. I freeze this in 1 ml aliquots and thaw and refreeze as needed.)
4. Perkin Elmer Amplitaq Gold
5. genomic DNA samples diluted to 50 ng/ μ l with sterile water

PCR reaction mixture:

15.8 μ l sterile water
2.5 μ l 10X PCR buffer
4 μ l dNTP premix
0.75 μ l primer #1 (Pdx1.5b or Pdx1.S5)
0.75 μ l primer #2 (Pdx1.S3 or Pdx1.3c)
1 μ l dil. DNA template
0.2 μ l Amplitaq Gold
25 μ l total volume

Cycling conditions:

1 cycle - 94°C x 6 min.

40 cycles - 94°C x 1 min., 60°C x 30 sec., 72°C x 45 sec.

1 cycle - 72°C x 7 min.

hold at 4°C.

Analysis of PCR products:

Load 10 μ l aliquots of reactions + 2 μ l of 6X gel loading buffer in a 1.5% mini-agarose gel (using a 15-well comb). Run gel approximately half-way down.

Example 3' reaction:

