

## PCR genotyping protocol for *Neurog3*<sup>HA.LCA</sup> mouse line

### PCR reagents:

1. oligo primers diluted to a 20 uM working concentration  
Ngn3HA-F 5'-TGCAGTGACCTCTAAGTCAGAGGCT-3'  
Ngn3HA-R 5'- ACCCACTTCTGCTTCGGAGCAGT-3'

Ngn3HA-F+ Ngn3HA-R yields a **230bp WT** (wild type allele) band and **330bp TM** (targeted mutant allele) band.

2. Perkin Elmer PCR buffer with MgCl<sub>2</sub>
3. 1.25mM dNTP premix (dNTP premix is made by using 100 mM NEB dNTP's. The premix contains 250 µl of each dNTP - A, C, G, &T - and 19 ml sterile water and is stable at -20°C. I freeze this in 1 ml aliquots and thaw and refreeze as needed.)
4. Perkin Elmer Amplitaq Gold
5. genomic DNA samples diluted to 50 ng/ul with Promega Nuclease-Free water

### PCR reaction mixture:

15.8 ul sterile water  
2.5 ul 10X PCR buffer  
4 ul dNTP premix  
0.75 ul primer Ngn3HA F  
0.75 ul primer Ngn3HA R  
1 ul dil. DNA template  
0.2 ul Amplitaq Gold  
25 ul total volume

### Cycling conditions:

1 cycle - 94° x 5 min.  
35 cycles - 95° x 30 sec., 59° x 30 sec., 72° x 30 sec.  
1 cycle - 72° x 7 min.  
hold at 4°C.

### Analysis of PCR products:

Load 10 ul aliquots of reactions + 2 ul of 6X gel loading buffer in a 1.5% mini-agarose gel (using a 15-well comb). Run gel approximately half-way down with a DNA size marker appropriate for distinguishing the amplicons you expect.

