

PCR protocol for genotyping *Gck*^{K414E} mouse line

PCR reagents:

1. oligo primers at 20 μ M

GK primers

GKP2 5'-TGT CTC AAT TTG CTG TGT CCT CCA-3' (top)

GKP8 5'-ATG TGT GAG TGT GCC AAT AT GAG T-3' (bottom)

(GKP2 + GKP8 yields a **636bp WT** (wild type allele) band and a **741bp TM** (targeted mutant allele) band.

2. Perkin Elmer PCR buffer with $MgCl_2$
3. 1.25mM dNTP premix (dNTP premix is made by using 100 mM NEB dNTP's. The premix contains 250 μ l of each dNTP - A, C, G, & T - and 19 ml sterile water and is stable at -20°C. I freeze this in 1 ml aliquots and thaw and refreeze as needed).
4. Perkin Elmer Amplitaq Gold
5. genomic DNA samples diluted to 50 ng/ μ l with sterile water

PCR reaction mixture:

15.5 μ l sterile water
2.5 μ l 10X PCR buffer
4 μ l dNTP premix
0.75 μ l primer GKP2
0.75 μ l primer GKP8
1 μ l dil. DNA template
0.2 μ l Amplitaq Gold
25 μ l total volume

Cycling conditions:

1 cycle - 94°C x 6 min.
40 cycles - 94°C x 1 min., 60°C x 30 sec., 72°C x 30 sec.
1 cycle - 72°C x 7 min.
hold at 4°C.

Analysis of PCR products:

Load 10 μ l aliquots of reactions + 2 μ l gel loading buffer in a 1% mini-agarose gel (using a 15-well comb). Run gel approximately half-way down.

