

1.2. loxP-2 PCR

Primers for loxP-2 PCR:

F2: 5'-CAATCCTATCATTTACGCCCTGC-3'

R2: 5'-GTTTCCAGTAGAAAGACAGGTGGT-3'

Expected PCR Product:

Wildtype: 227 bp

Targeted: 261 bp

Reaction Mix:


Component	x1
Mouse genomic DNA	1.5 µl
Forward primer (10 µM)	1.0 µl
Reverse primer (10 µM)	1.0 µl
Premix Taq Polymerase	12.5 µl
ddH ₂ O	9.0 µl
Total	25.0 µl

Cycling Condition:

Step	Temp.	Time	Cycles
Initial denaturation	94 °C	3 min	
Denaturation	94 °C	30 s	
Annealing	62 °C	35 s	35 x
Extension	72 °C	35 s	
Additional extension	72 °C	5 min	

Result

Seven pups (31#, 32#, 36#, 38#, 42#, 43# and 44#) from clone 2B5 were identified positive by PCR screening for loxP-2.



cyagen
We help you discover life

1.3. Neo-del PCR

Primers for Neo-del PCR:

F1: 5'-GCGCCATAACTTCGTATAGCAT-3'

R1: 5'-TTGGCTTCTTCTACTGGAGCTGTC-3'

Expected PCR Product:

Wildtype: N.A.

Targeted: 517 bp

Reaction Mix:


Component	x1
Mouse genomic DNA	1.5 µl
Forward primer (10 µM)	1.0 µl
Reverse primer (10 µM)	1.0 µl
Premix Taq Polymerase	12.5 µl
ddH ₂ O	9.0 µl
Total	25.0 µl

Cycling Condition:

Step	Temp.	Time	Cycles
Initial denaturation	94 °C	3 min	
Denaturation	94 °C	30 s	
Annealing	62 °C	35 s	35 x
Extension	72 °C	35 s	
Additional extension	72 °C	5 min	

Result

Seven pups (31#, 32#, 36#, 38#, 42#, 43# and 44#) from clone 2B5 were identified positive by PCR screening for Neo-del.




1.4. PCR Result:

Seven pups (31#, 32#, 36#, 38#, 42#, 43# and 44#) from clone 2B5 were identified positive by PCR screening for loxP-2 and Neo-del, the positive pups were reconfirmed by PCR screening for Neo-del.



1.5. Suggested Breeding and Genotyping Assay for Tissue-specific Knockout Mice Generation


Wildtype allele




Mutant allele 1 (Targeted allele)




Mutant allele 2 (After Neo deletion)



Mutant allele 3 (After Cre recombination)



Legends



Step 1: Inter-cross heterozygous targeted mice to generate homozygous targeted mice

Primers for targeted allele:

F2: 5'-CAATCCTATCATTTACGCCCTGC-3'

R2: 5'-GTTTCCAGTAGAAAGACAGGTGGT-3'

Wildtype: 227 bp

Homozygotes: 261 bp

Heterozygotes: 261 bp/227 bp

Step 2: Breed a homozygous targeted mouse with a tissue-specific Cre delete mouse to generate mice that are heterozygous for a targeted allele and a hemizygous/heterozygous for the Cre transgene

Primers for targeted allele:

F2: 5'-CAATCCTATCATTTACGCCCTGC-3'

R2: 5'-GTTTCCAGTAGAAAGACAGGTGGT-3'

Heterozygotes: 261 bp/227 bp

Primers for Cre transgene:

Forward1: 5'-CATATTGGCAGAACGAAAACGC-3'

Reverse1: 5'-CCTGTTTCACTATCCAGGTTACGG-3'

Cre amplicon: 413 bp

Step 3: Breed heterozygous, Cre⁺ mice with homozygous mice. Approximately 25% of the progeny from this mating will be homozygous for the targeted allele and hemizygous/heterozygous for the Cre transgene. The pups can be screened by the same assay as described above. The tissue-specific gene deletion can be confirmed by the following primers:

Primers for targeted allele:

F2: 5'-CAATCCTATCATTTACGCCCTGC-3'

R2: 5'-GTTTCCAGTAGAAAGACAGGTGGT-3'

Conditional KO allele: 261 bp

Wildtype allele: 227 bp

F3: 5'-TGCTTGGCAGCAAACCTATAAATGG-3'

R3: 5'-TTCCATTACCCTAGAACTGGCTTC-3'

Constitutive KO allele: 456 bp