



# Vanderbilt Genome Editing Resource

## GENOTYPING PROTOCOL: *SOX17<sup>em3Mgn</sup>*

**Investigator:** Mark Magnuson, VICC and DRTC

**Genome Edit:** Sox17 CR2 -50 bp deletion

**Allele Name:** *Sox17<sup>em3Mgn</sup>*

**Primers:**

Sox17CR2Fwd1: CTTGGACTTGTTCTTCAATCTTCC

Sox17CR2Rev2: CCAGAGGAACTCGTAAAGCTG

**PCR products:**

WT = 394 bp

Sox17 -50 bp allele = 344 bp

Components	50 ul reaction	Final concentration	PCR program*
5 x Phusion Reaction Buffer (NEB #M0530S)	10.0 µL	1X	98°C, 30 seconds
10 mM dNTPs	1.0 µL	200 µM	98°C, 10 seconds
10 µM Sox17CR2Fwd1	2.5 µL	0.5 µM	62°C, 10 seconds
10 µM Sox17CR2Rev2	2.5 µL	0.5 µM	72°C, 40 seconds
Phusion HF DNA Polymerase	0.5 µL	0.02 U/µl	Go to 2, 35 X
Nuclease-free water	32.5 µL		72°C, 2 minutes
Genomic DNA	1.0 µL	Less than 1 µg, variable	4°C, ∞

